

FOR EMERGENCY USE!



Cochlear Implant Users need
accessibility to communication
and Cochlear Implant Care & help
during hospitalization

NAME: _____

DATE OF BIRTH: _____

ADDRESS: _____

MY CONTACT FOR QUESTIONS: _____

I HAVE HEARING LOSS AND WEAR THESE DEVICES:

Left: Cochlear Implant / Hearing Aid

Model and Brand: _____

Right: Cochlear Implant / Hearing Aid

Model and Brand: _____

Important: Do not place devices on the wrong ears!



Cochlear implants (CI) replace natural hearing with electronic hearing. Keep in mind:

I cannot hear anything without my external Cochlear Implant (CI) processor(s).


Look at me! Eye contact and the possibility to lip-read will help us to communicate.

CI processor on and off: Attach/remove batteries. Alternatively, push the lower button briefly to turn on and hold for 2 seconds to turn off.

To attach the processor: Connect the coil of the external processor to the invisible magnet under the scalp (located slightly above / behind the ear).

To replace rechargeable batteries: Turn and remove battery from the processor. Put it on the charger until the LED is green.

To replace disposable batteries: Slide the battery cover off to expose the disposable batteries. **Batteries usually last:** _____



At night, place the cochlear implant processor(s) in the dehumidifier if provided.

Please do not touch the remote control(s).

DO NOT perform MRI without prior clarification and professional agreement with the otolaryngologist at the cochlear implantation clinic and the **manufacturer**:

If necessary, please write down important messages clearly so that I can read them!

© EURO-CIU, asbl www.eurociu.eu

© Bayerischer Cochlea-Implantat-Verband e.V., www.bayciv.de