

FOR EMERGENCY USE!

NAME: _____

DATE OF BIRTH: _____

ADDRESS: _____

MY CONTACT FOR QUESTIONS: _____



Cochlear Implant Users need
accessibility to communication
and Cochlear Implant Care & help
during hospitalization

I HAVE HEARING LOSS AND WEAR THESE DEVICES:

Left: Cochlear Implant / Hearing Aid

Model and Brand: _____

Right: Cochlear Implant / Hearing Aid

Model and Brand: _____

Important: Do not place devices on the wrong ears!



Cochlear implants (CI) replace natural hearing with electronic hearing. Keep in mind:

I cannot hear anything without my external Cochlear Implant (CI) processor(s).

Look at me! Eye contact and the possibility to lip-read will help us to communicate.

CI processor on and off: Attach or remove batteries.

To attach the processor: Connect the coil of the external processor to the invisible magnet under the scalp (located slightly above / behind the ear).



To replace rechargeable batteries: Rotate and withdraw the battery from the processor and attach it to the charger until the LED turns green.

To replace disposable batteries: Slide the battery cover downwards to expose the disposable batteries. **Batteries usually last:** _____

At night, place the cochlear implant processor(s) in the dehumidifier if provided.

Please do not touch the remote control(s).

DO NOT perform MRI without prior clarification and professional agreement with the otolaryngologist at the cochlear implantation clinic and the **manufacturer**:

If necessary, please write down important messages clearly so that I can read them!

© EURO-CIU, asbl www.eurociu.eu

© Bayerischer Cochlea-Implantat-Verband e.V., www.bayciv.de