

FOR EMERGENCY USE!

NAME: _____

DATE OF BIRTH: _____

ADDRESS: _____

MY CONTACT FOR QUESTIONS: _____



Cochlear Implant Users need
accessibility to communication
and Cochlear Implant Care & help
during hospitalization

I HAVE HEARING LOSS AND WEAR THESE DEVICES:


Left: Cochlear Implant / Hearing Aid

Model and Brand: _____

Right: Cochlear Implant / Hearing Aid

Model and Brand: _____

Important: Do not place devices on the wrong ears!



Cochlear implants (CI) replace natural hearing with electronic hearing. Keep in mind:

I cannot hear anything without my external Cochlear Implant (CI) processor(s).

Look at me! Eye contact and the possibility to lip-read will help us to communicate.

CI processor on and off: Push the red button on the Rondo2.

To attach the Rondo2 processor: Connect the external processor to the invisible magnet under the scalp (located slightly above / behind the ear).

To charge batteries: Place the Rondo2 onto the centre of the charging station with the MED-EL logo pointing upwards.

Batteries usually last: _____



Please do not touch the remote control(s).

DO NOT perform MRI without prior clarification and professional agreement with the otolaryngologist at the cochlear implantation clinic and the **manufacturer**:

If necessary, please write down important messages clearly so that I can read them!

© EURO-CIU, asbl www.eurociu.eu

© Bayerischer Cochlea-Implantat-Verband e.V., www.bayciv.de